

Development of iPSC-derived Alzheimer Disease (AD) platforms for phenotypic screening

Angela C. Murchison¹, Peng Zhou¹, Christopher Noel¹, Daniel Haag¹, Ji W¹, Tao Huang¹, and Wayne W. Poon^{1,2}


¹NeuCyte, Inc., 319 N. Bernardo Avenue, Mountain View CA 94043

²Institute for Memory Impairments and Neurological Disorders, UC Irvine, Irvine, CA 92697


Abstract

Here we describe the development of iPSC-derived platforms of AD that can be used to screen for preclinical targets and drug candidates. The foundation of these platforms utilizes a proprietary method to generate pure populations of either NGN2-excitatory or ASCL1/DLX2-inhibitory neurons. We developed a neuron/astrocyte co-culture platform to uncover *APOE4*-specific electrophysiological phenotypes in which *APOE4* is the greatest genetic risk factor of sporadic AD. This platform can be complemented with isogenic iPSC-derived astrocytes and microglia to interrogate both cell autonomous and non-cell autonomous mechanisms for genetic risk factors that either confer increased or decreased AD risk. We also report the development of an iPSC-based transgenic tau overexpression platform to investigate how different tau isoforms promote prion-like tau spreading and be used to screen for compounds to mitigate tau spreading. Tau overexpression causes neuron-neuron tau transmission and leads to the accumulation of key neuropathological tau biomarkers. Taken together, these innovative phenotypic assays can be translated to high-throughput drug screening efforts for AD.


Induced Pluripotent Stem Cells Facilitate Translational Disease Modeling for Drug Discovery





NeuCyte SynFire® Platform Enables Development of Phenotypic Assays for Drug Screening



APOE Isogenic iPSCs Facilitate Human-specific Phenotypic Discovery




Development of an Electrophysiological MEA Phenotypic Assay for *APOE4*



(A) Representative immunocytochemical images of induced excitatory and inhibitory *APOE4* and *APOE3* neurons cultured on primary astrocytes (DPP14). **(B)** Representative raster plots demonstrating the normal development of mature neuronal networks over 21 days. **(C)** *APOE4* excitatory neurons exhibit impaired network synchrony.

Development of an iPSC-based Tau Pathological Screening Assay for Drug Discovery



Future Directions

- Incorporate iPSC-derived Astrocytes and Microglia into Co-cultures for investigating Neuroinflammatory mechanisms
- Develop 3D organoid/assembleid screening assays
- Translate Phenotypic assays to Drug Screening and Drug Discovery
- Examine *APOE4*/tau/GABA tripartite relationship