

Human iPSC-derived sensory neurons and their potential as an alternative to primary DRG cells for drug discovery



Minami Hiranuma¹, Yuichi Okuda¹, Tomohisa Watanabe¹, and Robert Annand²

1. REPROCELL, Yokohama, Japan
2. REPROCELL, Beltsville, MD USA

www.reprocell.com

Sensory nerves transmit pain signals generated by trauma, temperature, or chemicals to the central nervous system. Primary afferent nerves, evoked by non-nociceptive and nociceptive stimuli, reside in the dorsal root ganglion (DRG) of the spinal cord. Primary DRG cells have been used as an *in vitro* model of nociceptive responses for decades. However, they are low-yielding, challenging to culture, and difficult to access. To negate accessibility issues and improve ease of culture, REPROCELL has developed an alternative sensory

neuron model using human-derived induced pluripotent stem cells (iPSC). This model was analyzed to confirm its phenotypical relevance and functionality. Immunocytochemistry and real-time PCR were used to compare the expression of sensory markers to primary DRG cells. Multi-Electrode Array (MEA) analysis was used to confirm whether these induced hiPSC-derived sensory nerves display the same functionality as primary DRG cells.

Conclusion: The sensory nerves we have developed in this study express sensory neuron-related proteins and genes, and responded to capsaicin, menthol, temperature, and bradykinin. We have therefore concluded that these iPSC-derived sensory neurons can replace primary DRG cells as a model for drug discovery research.

1 Expression of sensory neuron markers




Fig 1. Immunofluorescence staining day 14 after thawing and seeding. Neurons were seeded at 1.5×10^5 cells/well (24 well plates) in Sensory Neuron Culture Medium. They were stained for neural markers, peripheral nerve markers, and various sensory nerve-related proteins including: Tuj1 (mature nerve marker), Peripherin (peripheral nerve marker), Brn3a (sensory nerve marker), TRPV1 (capsaicin receptor), TRPM8 (menthol receptor), Nav 1.7 (Na channel), TRKA (NGF receptor), TRKB (BDNF receptor), and TRKC (neurotrophin-3 receptor).

4-8 Electrophysiological analysis by MEA.





Fig 4-8. Electrophysiological analysis by MEA (n=3 wells). Electrophysiological response of sensory nerves to (4) 100 nM capsaicin 27 days after seeding, (5) 1nM menthol 32 days after seeding, (6) 100 nM Bradykinin 39 days after seeding, (7) temperature 28 days after seeding, (8) electrical stimulation. For Fig 7, temperature was increased from 37°C to 46°C at 28 days after seeding. For Fig 8, pulse conditions used were 1.2 V, 750 µs/phase, 40 times, 300 Hz. Cells seeded in CytoView MEA 24-White at 1.5×10^5 cells/well/1.1 mm². Maestro Edge (Axion) used for MEA analysis. (A) Raster Plot, (B) Mean Firing Rate (Hz), (C) Number of Bursts.

2 Creation of Sensory Neurons



9 Electrophysiological response of sensory nerves to ProTX-II




Fig 9. Electrophysiological response of sensory nerves to 1 µM ProTX-II (n=3 wells). Cells seeded in CytoView MEA 24-White at 1.5×10^5 cells/well/1.1 mm². Maestro Edge (Axion) used for MEA analysis. (A) Raster plot before addition, (B) Raster plot 30 min after addition when exposed to 1 µM ProTX-II at 43 days after thaw seeding, (C) Mean Firing Rate (Hz), (D) Number of Bursts.

3 Expression of sensory neuron genes (n= 1)




Fig 3. Real-time PCR analysis on day 14 after thawing and seeding. Gene expression of Peripherin, Brn3a, TRP channel, and high-affinity neurotrophin receptors was increased significantly in DRG compared with hiPSCs. There was no marked difference in Brn3a, TRPV1, TRPM8, and TRKB gene expression.