

Instruction manual

ReproMed™ iPSC Medium

Cat. No. RCRM101

Version 1.3

Overview

This instruction manual explains how to passage and maintain human iPSC which have been cultured on feeder cells to feeder-free culture using **ReproMed™ iPSC Medium** (after supplementing with bFGF).

Please use this product in a biosafety cabinet.

About product

Sale of this product to a third party or other commercial purposes is strictly prohibited without prior permission from ReproCELL.

Storage

This product is shipped frozen. Upon receipt store immediately at $-20\text{ }^{\circ}\text{C}$. Thaw overnight at $2\text{ }^{\circ}\text{C}$ to $8\text{ }^{\circ}\text{C}$ before use. Aliquot and store at $-20\text{ }^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Before using, supplement with bFGF and use within 2 weeks.

Characteristics

- Feeder-free culture.
- Each lot has been tested on PChiPS771 human iPSC established by ReproCELL.
- Quality control testing of critical criteria, including osmolality, pH, sterility, and mycoplasma has been performed on each lot.
- Does not contain any serum.

Please take note of the following points when using this product.

- It is recommended to incubate cells in a **5% O₂, 5% CO₂ incubator**.
(It is possible to use a 5% CO₂ incubator, but the proliferation rate is about twice as fast in a 5% O₂, 5% CO₂ incubator.)
- Supplement with bFGF (at a final concentration of 10 ng/mL) before using. Hereinafter, this product is referred to as **ReproMed™ iPSC Medium**.
- Allow ReproMed™ iPSC Medium to sit for 30 min at room temperature to allow the medium to warm up to room temperature before using. Do not use a water bath to warm the medium.
- While 10 ng/mL of bFGF is recommended, the amount of bFGF can be adjusted to suit different cells.
- After the initial transition to feeder-free culture, there may be many floating cells observed in the medium. This is normal and will decrease with subsequent passages.
- As listed below are the recommended dissociation solutions and coating reagent for feeder-free culture:
Dissociation solutions: 0.5mM EDTA/PBS+ ESGRO Complete™ Accutase™ & TrypLE™ Select Enzyme (1X), no phenol red
Coating reagent: iMatrix-511

Table of contents

- **Required reagents & equipment**
- **ReproMed™ iPSC Medium culture protocol**
 - **1. Coating**
 - **2. Transition from feeder-dependent to feeder-free culture using ReproMed™ iPSC Medium**
 - **3. ReproMed™ iPSC Medium culture protocol**
 - **4. Morphology of feeder-free human iPSC grown in ReproMed™ iPSC Medium**
- **Frequently asked questions**

Product

Product Description	Cat. No.	Amount	Storage
ReproMed™ iPSC Medium	RCRM101	500 mL	-20°C

Required reagents & equipment

Product Description	Cat. No.	Amount	Storage
bFGF	REPROCELL RCHEOT002	25 µg	-80°C
iMatrix-511	REPROCELL NP892-011	350 µg	4°C
Dissociation Solution for human ES/iPSC	RCHETP002	30 mL	-20°C
TrypLE™ Select Enzyme (1X), no phenol red	GIBCO Cat. 12563011	100 mL	room temperature
ESGRO Complete™ Accutase™	Millipore SF006	100 mL	4 °C
UltraPure™ 0.5M EDTA, pH 8.0	Invitrogen 15575020	100 mL	room temperature
Y27632	Wako 257-00511	-	-
PBS (-). Ca ²⁺ - and Mg ²⁺ -Free	-	-	-
60 mm Tissue Culture Dish	-	-	-
5% CO ₂ incubator	-	-	-

ReproMed™ iPSC Medium culture protocol

1. Coating

<For 4 x 60 mm dish (surface area = 22 cm²)>

- Place iMatrix-511 on ice.
- Aliquot 8 mL PBS (-) to a 15 mL tube and place it on ice.
- Aliquot 88 µL of iMatrix-511 to the PBS (-) and mix gently. (92-fold dilution of iMatrix-511)
- Coat each dish with 2 mL of the prepared iMatrix-511 solution. (0.5 µg iMatrix-511/cm²)
- Incubate for 1 hour in a 37 °C, 5% CO₂ incubator. Seal the coated dish with parafilm and store at 4 °C if not used immediately and use within a week.

2. Transition from feeder-dependent culture on 100 mm dish to feeder-free culture using ReproMed™ iPSC Medium on 60 mm dish

Note 1: The reagent volumes stated below are for transition to a 60 mm dish.

Note 2: Use a 100 mm dish that is about 30% confluent for the transition (Refer to Figure 1).




Figure 1
Alkaline Phosphatase (ALP) stained 30% confluent 100mm dish

- Aspirate to remove the spent medium.
 - Wash dish with 4 mL of PBS (-).
 - Add 1.5 mL of dissociation solution for human ES/iPSC.
 - Incubate for 4 min in a 37 °C, 5% CO₂ incubator.
 - Aspirate to remove the dissociation solution for human ES/iPS cells.
 - Wash dish with 2 mL of PBS (-).
 - Add 1.5 mL of 0.5X TrypLE Select (TrypLE™ Select Enzyme (1X) diluted with an equal volume of 0.5 mM EDTA).
 - Incubate for 7 min in a 37 °C, 5% CO₂ incubator.
 - Flush the surface of the dish 10 times with a P1000 pipette to detach the cells from the dish.
 - Collect and transfer the cells to a 15 mL tube.
 - Collect the remaining cells with 2 mL/dish of ReproMed™ iPSC Medium containing Y27632 (at a final concentration of 10 µM).
- Note 3:** Adding Y27632 increases cell viability and stability of the culture.
- Centrifuge at 300×g for 5min.
 - Aspirate to remove the supernatant.
 - Resuspend the cells with 2 mL of ReproMed™ iPSC Medium with 10 µM Y27632 and count the

number of viable cells.

15. Aliquot an appropriate volume of cells into a new 15 mL tube and create a 4 mL cell suspension with ReproMed™ iPSC Medium with 10 μM Y27632. Transfer this cell suspension to an iMatrix-511 coated 60 mm dish.

(※) Recommended seeding number:

For the first 3 passages after the transition to feeder-free culture, seed the cells at 5×10^5 cells/dish.

After 3 passages, optimize according to cell type.

Option (1): For slow growing strains:

5×10^5 cells/dish

Option (2): For fast growing strains with a high ability to remain undifferentiated:

2×10^5 cells/dish

16. Incubate in a 37 °C, 5% CO₂ incubator.
17. Replace with 4 mL/dish of fresh ReproMed™ iPSC Medium (without Y27632) daily.

Note 4: Option (1): Replace with 4 mL/dish of fresh ReproMed™ iPSC Medium (without Y27632) on Day 1, 24 hours after seeding.

Option (2): Replace with 4 mL/dish of fresh ReproMed™ iPSC Medium with 10 μM Y27632 on Day 1, 24 hours after seeding. On Day 2, replace with 4 mL/dish of fresh ReproMed™ iPSC Medium (without Y27632), 48 hours after seeding. (Including Y27632 in the culture medium for 24 hours increases the stability of the cells cultured.)

Note 5: Passage the cells once they have reached 80~90% confluent. (Cells seeded on Day 0 can be passaged on Day 7. For cells that are grown in a 5% O₂, 5% CO₂ incubator, cells can be passaged on Day 5.)

Note 6: After the initial transition to feeder-free culture, there may be many floating cells observed in the medium. This is normal and will decrease with subsequent passages.

3. ReproMed™ iPSC Medium culture protocol

Note 7: The reagent volumes stated below are for cultures in a 60 mm dish.

Note 8: The protocol below is different from the transition of feeder-dependent to feeder-free culture protocol.

1. Aspirate to remove the spent ReproMed™ iPSC Medium.
2. Wash twice with 2 mL of PBS (-).
3. Add 1 mL/ dish of 0.5 mM EDTA/PBS and incubate for 2 min at room temperature on a rotator before removing the 0.5 mM EDTA/PBS.
4. Add 1 mL/dish of ESGRO Complete™ Accutase™.
5. Incubate for 7 min in a 37 °C, 5% CO₂ incubator.
6. Flush the surface of the dish with a P1000 pipette to detach the cells from the dish.

7. Collect and transfer the cells to a 15 mL tube.
8. Collect the remaining cells with 1 mL/dish of ReproMed™ iPSC Medium with 10 μM Y27632.

Note 9: Adding Y27632 increases cell viability and stability of the culture.

9. Centrifuge at 300×g for 5 min.
10. Aspirate to remove the supernatant.
11. Resuspend the cells with 2 mL of ReproMed™ iPSC Medium with 10 μM Y27632 and count the number of viable cells.

12. Aliquot an appropriate volume of cells into a new 15 mL tube and create a 4 mL cell suspension with ReproMed™ iPSC Medium with 10 μM Y27632. Transfer this cell suspension to an iMatrix-511 coated 60 mm dish.

(※) Recommended seeding number:

For the first 3 passages after the transition to feeder-free culture, seed the cells at 5×10^5 cells/dish.

After 3 passages, optimize according to cell type.

Option (1): For slow growing strains:

5×10^5 cells/dish

Option (2): For fast growing strains with a high ability to remain undifferentiated:

2×10^5 cells/dish

Incubate in a 37 °C, 5% CO₂ incubator.


13. Replace with 4 mL/dish of fresh ReproMed™ iPSC Medium (without Y27632) daily.

Note 10: Option (1): Replace with 4 mL/dish of fresh ReproMed™ iPSC Medium (without Y27632) on Day 1, 24 hours after seeding.

Option (2): Replace with 4 mL/dish of fresh ReproMed™ iPSC Medium with 10 μM Y27632 on Day 1, 24 hours after seeding. On Day 2, replace with 4 mL/dish of fresh ReproMed™ iPSC Medium (without Y27632), 48 hours after seeding. (Including Y27632 in the culture medium for 24 hours increases the stability of the cells cultured.)

Note 11: Passage the cells once they have reached 80~90% confluent. (Cells seeded on Day 0 can be passaged on Day 7. For cells that are grown in a 5% O₂, 5% CO₂ incubator, cells can be passaged on Day 5.)

4. Morphology of feeder-free human iPSC (253G1) grown in ReproMed™ iPSC Medium (seeded at 5×10^5 cells/dish and cultured in a 5% CO₂ incubator)



Frequently asked questions

Q1. How do I switch to ReproMed™ iPSC Medium from another feeder-free medium?

A1. Switch to ReproMed™ iPSC Medium when passaging cells. After dissociating cells, centrifuge cells and aspirate to remove the supernatant. Resuspend the cells in ReproMed™ iPSC Medium for counting and seeding.

REPROCELL Inc.

<https://www.reprocell.com>

E-mail: info_en@reprocell.com